

BPX5

A Capillary Column which will perform 80% of all applications

What is BPX5?

BPX5 is a 5% phenyl (equivalent) / 95% methyl polysilphenylene / siloxane phase.

Why use a 5% phenyl phase?

It is estimated that 80% of all analyses can be carried out on one column with one set of dimensions and film thickness. The 5% phenyl column is still a non-polar column, but the small amount of phenyl content gives the phase a little edge in selectivity compared with the 100% methyl phases.

The non-polar phases tend to have lower bleed and are more robust compared with polar phases.

What can the BPX5 column be used for?

The BPX5 is a multi-purpose column which can be used for a range of applications from non-polar components (such as those in gasoline) to polar components such as amines, pesticides, triglycerides, drugs and polycyclic aromatic hydrocarbons.

Equivalent Columns

If you are using any 5% phenyl column, BPX5 can replace it. The table below is an excellent replacement reference.

This column is a direct replacement for:

CP-SIL 8CB	OV-5ms
CP-SIL 8CBms	PE-5
DB-5	PTE-5
DB-5.625	Rtx-5
DB-5-MS	Rtx-5MS
HP-5	SPB-5
HP-5MS	VB-5
HP-5TA	XTI-5
HP-Ultra2	ZB-5 (RH-5)
MDN-5, MDN-5S	007-5
OV-5	007-5MS

Operating Temperatures


0.25µm – 1.0µm -40°C to 360°C/370°C
>1.0µm -40°C to 350°C/360°C

PHENOLS

USEPA 625 Phenols Mix

Phase: BPX5, 0.25mm
Column: 30m x 0.25mm ID
Injector Temp.: 240°C
Injector Mode: Split, 80:1, 1.3ng O.C.
Injection Volume: 0.5mL
Carrier Gas: Helium, 25cm/sec
Initial Oven Temp.: 40°C, 1min
Rate 1: 10°C/min
Temp. 1: 140°C, 0min
Rate 2: 20°C/min
Temp.: 280°C, 0min
Detector: FID, 300°C
Part No.: 054101


- | | |
|----------------------------|--------------------------------|
| 1. Phenol | 7. 2,4,6-Trichlorophenol |
| 2. 2-Chlorophenol | 8. 2,4-Dinitrophenol |
| 3. 2-Nitrophenol | 9. 4-Nitrophenol |
| 4. 2,4-Dimethylphenol | 10. 2-Methyl-4,6-dinitrophenol |
| 5. 2,4-Dichlorophenol | 11. Pentachlorophenol |
| 6. 4-Chloro-3-methylphenol | |



PESTICIDES - BROAD RANGE PESTICIDE MIXTURE

Phase: BPX5, 0.25µm
Column: 30m x 0.22mm
 Initial Temp.: 50°C, 1min
 Rate 1: 30°C/min
 Temp. 2: 190°C
 Rate 2: 4°C/min
 Final Temp.: 280°C, 5min
 Detector: 5973 MSD
 Injector Mode: Splitless, 1.0min
 Carrier Gas: Helium, 35cm/sec
 Injection Volume: 1µL
 Interface/Source/Quad: 320°C, 250°C, 150°C
Part No.: 054142

- | | |
|-----------------------------|---------------------|
| 1. Methyl-Dymron decomp. I | 15. Chlorpyrifos |
| 2. Methyl-Dymron decomp. II | 16. Pendimethalin |
| 3. Trichlorfon (DEP) | 17. Isofenphos |
| 4. Etridiazol | 18. Methyl-Dymron |
| 5. Chloroneb | 19. Captan |
| 6. Benfluralin | 20. Butamifos |
| 7. Pencycuron | 21. Napropamide |
| 8. Simazine (CAT) | 22. Flutolanil |
| 9. Diazinon | 23. Isoprothioate |
| 10. Propyzamide | 24. Isoxathion |
| 11. Chlorothalonil (TPN) | 25. Bensulide |
| 12. Terbutcarb (MBPMC) | 26. Mepronil |
| 13. Tolcolphos-Methyl | 27. Pyridaphenthion |
| 14. Fenitrothion (MEP) | 28. Iprodione |




Note: The non-polar nature of the BPX5 and the high maximum operating temperature are the main reasons our customers prefer this column.

ANALYSIS OF CHLORINE, PHOSPHOROUS AND NITROGEN BASED PESTICIDES

Phase: BPX5, 0.25µm
Column: 30m x 0.25mm ID
 Initial Temp.: 50°C, 1min
 Rate 1: 4°C/min
 Final Temp.: 280°C, 1.5min
 Detector: 5973 MSD
 Injector Mode: Splitless, 0.50min
 Carrier Gas: Helium, 33cm/sec
 Injection Volume: 0.5µL
Part No.: 054101

- | | |
|-----------------------|---------------------|
| 1. DDVP | 25. Thiobancarb |
| 2. Etridiazol | 26. Fthalide |
| 3. Chloroneb | 27. Pendimethalin |
| 4. Molinate | 28. Isofenphos |
| 5. BPMC | 29. Methyl-Dymron |
| 6. Benfluralin | 30. Captan |
| 7. Pencycuron | 31. Butamifos |
| 8. Simazin | 32. Napropamide |
| 9. Diazinon-ox | 33. Flutolanil |
| 10. Propyzamide | 34. Prethirachlor |
| 11. Diazinon | 35. Isoxathon-ox |
| 12. TPN | 36. Isoprothioate |
| 13. IBP | 37. Tricyclazole |
| 14. ECP | 38. Buprofezine |
| 15. MBPMC | 39. Isoxathion |
| 16. Bromobutide | 40. Mepron:1 |
| 17. MEP-ox | 41. Chloronitrofen |
| 18. Tolclophos-Methyl | 42. EDDP |
| 19. Simetryn | 43. EPN-ox |
| 20. MEP | 44. Pyridaphenthion |
| 21. Esprocarb | 45. Iprodione |
| 22. Malathion | 46. EPN |
| 23. Probenazole | 47. Bifenox |
| 24. Chlorpyrifos | 48. Mefenacet |




Note: The BPX5 column is an ideal partner column to BP10 (1701 phase), BPX608, or BPX50.

PESTICIDES - EPA METHOD 608 AND ADDITIONAL PESTICIDES

Phase: BPX5, 0.25µm
Column: 30m x 0.22mm ID
 Initial Temp.: 160°C
 Rate 1: 7°C/min
 Final Temp.: 290°C, 5min
 Detector: ECD, 320°C
 Injector Mode: Split, 300:1
 Carrier Gas: Helium, 25cm/sec
 Injection Volume: 0.5µL
Part No.: 054142


- | | |
|------------------------|------------------------|
| 1. Etridiazol | 14. gamma-Chlordane |
| 2. Chloroneb | 15. alpha-Chlordane |
| 3. Trifluralin | 16. Endosulfan I |
| 4. alpha-BHC | 17. Dieldrin |
| 5. gamma-BHC | 18. p,p'-DDE |
| 6. beta-BHC | 19. Endrin |
| 7. Chlorothalonil | 20. Endosulfan II |
| 8. omega-BHC | 21. p,p'-DDD |
| 9. Heptachlor | 22. Endrin aldehyde |
| 10. Chlorpyrifos | 23. p,p'-DDT |
| 11. Aldrin | 24. Endosulfan sulfate |
| 12. Heptachlorepoixide | 25. Endrin ketone |
| 13. Captan | 26. Methoxychlor |
| | 27. Mirex |

Note: BPX5 column exhibits very low column bleed with selective detectors like ECD.




USEPA METHOD 608.1 - ORGANOCHLORINE PESTICIDES

Phase: BPX5, 0.25µm
Column: 30m x 0.22mm ID
 Initial Temp.: 50°C, 1min
 Rate 1: 30°C/min
 Temp. 2: 150°C
 Rate 2: 6°C/min
 Final Temp.: 260°C
 Detector: 5973 MSD
 Injector Mode: Splitless, 0.5min
 Carrier Gas: Helium, 30cm/sec
 Injection Volume: 1µL
 Interface/Source/Quad: 320°C, 250°C, 150°C
Part No.: 054142



USEPA METHOD 614 - ORGANOPHOSPHOROUS PESTICIDES


Phase: BPX5, 0.25µm
Column: 30m x 0.22mm ID
 Initial Temp.: 50°C, 1min
 Rate 1: 30°C/min
 Temp. 2: 150°C
 Rate 2: 10°C/min
 Final Temp.: 320°C
 Detector: 5973 MSD
 Injector Mode: Splitless, 0.5 min
 Carrier Gas: Helium, 30cm/sec
 Injection Volume: 1µL
 Interface/Source/Quad: 320°C, 250°C, 150°C
Part No.: 054142



TRIGLYCERIDES - ANALYSIS OF MONO-, DI-, AND TRIGLYCERIDES

Phase: BPX5, 0.25 μ m
Column: 12m x 0.32mm ID
Initial Temp.: 80°C, 1min
Rate: 10°C/min
Final Temp.: 365°C, 10min
Detector: FID
Injector Mode: On-Column (SGE OCI-5)
Carrier: He, 10psi
Part No.: 054118


Note: Triglycerides require high oven temperatures for elution. Therefore, the BPX5 column is the ideal choice for this analysis.



TRIGLYCERIDE DISTRIBUTION IN MILK FAT

Phase: BPX5, 0.25 μ m
Column: 12m x 0.53mm ID
Initial Temp.: 100°C, 0.5 min
Rate 1: 50°C/min
Temp.: 280°C
Rate 2: 10°C/min
Final Temp.: 360°C, 5 min
Carrier Gas: Helium, 6psi
Injection Mode: On-column (SGE OCI)
Injection Volume: 1 μ L
Part No.: 054133


Note: The low bleed nature of this column at high oven temperatures ensures reproducible analyses.



PEAK No.	TRIGLYCERIDE CARBON No.	COMPOSITION (%)
1	T24	0.12
2	T26	0.41
3	T28	0.64
4	T30	1.30
5	T32	2.61
6	T34	5.99
7	T36	11.22
8	T38	13.49
9	T40	10.55
10	T42	7.36
11	T44	6.90
12	T46	7.10
13	T48	8.41
14	T50	10.49
15	T52	9.34
16	T54	3.95
17	T56	0.22

STANDARD MIXTURE OF TRIGLYCERIDES

Phase: BPX5, 0.25 μ m
Column: 12m x 0.53mm ID
Initial Temp.: 100°C, 0.5min
Rate 1: 50°C/min
Temp.: 280°C
Rate 2: 10°C/min
Final Temp.: 360°C, 5min
Carrier Gas: Helium, 6psi
Detector: FID
Injector Mode: On-Column (SGE OCI-5)
Injection Volume: 1 μ L
Part No.: 054133




PEAK No.	TRIGLYCERIDE CARBON No.
1	T36
2	T38
3	T40
4	T42
5	T44
6	T46
7	T48
8	T50
9	T52
10	T54

AMINES - ANALYSIS OF ANILINE IMPURITIES

Phase: BPX5, 1.0 μ m
Column: 12m x 0.32mm ID
 Initial Temp.: 50°C, 1min
 Rate 1: 15°C/min
 Temp. 1: 120°C
 Rate 2: 20°C/min
 Final Temp.: 330°C, 10min
 Detector: FID, 380°
 Injector: Split
 Carrier Gas: He, 4psi
Part No.: 054127

Note: Thick film BPX5 column provides excellent separation with high temperature compatibility.




1. Methanol
2. Cyclohexane
3. Cyclohexanone
4. Cyclohexylamine
5. Cyclohexanol
6. Aniline
7. Toluidine
8. Unknown
9. Diaminodiphenylmethane
10. p,p'-Diaminodiphenyl-methane
11. Unknown

ANALYSIS OF AROMATIC AMINES

Phase: BPX5, 0.25mm
Column: 30m x 0.25mm ID
 Injector Temp.: 240°C
 Injector Mode: Split, 80:1, 3.25ng On Col
 Injection Volume: 0.5 μ L
 Carrier Gas: Helium, 30cm/sec @ 40°C
 Initial Oven Temp.: 40°C, 1min
 Rate 1: 20°C/min
 Final Temp.: 240°C, 0min
 Detector: FID, 300°C

Part No.: 054101

Note: Every BPX5 column is supplied with a test chromatogram which is produced from that column.




1. Aniline
2. 2-Methylpiperidine
3. N,N-Dimethylaniline
4. Benzylamine
5. α -Phenylethylamine
6. 2,6-Dimethylaniline
7. N,N-Diethylaniline

FRAGRANCE - ANALYSIS OF LAVENDER ESSENTIAL OIL

Phase: BPX5, 0.25 μ m
Column: 30m x 0.25mm ID
 Initial Oven Temp.: 45°C, 1 min
 Rate 1: 5°C/min
 Final Temp.: 250°C, 15 min
 Detector: 5973 MSD
 Injector Mode: Split (50:1)
 Injector Temp.: 250°C
 Carrier Gas: Helium, 35 cm/sec @ 45°C
 Injection Volume: 1 μ L
 Concentration (ng/mL): 1:20 dilution in DCM
 GC model: HP 6890 GC, 5973 MSD
 Interface/Source/Quad: 320°C, 250°C, 106°C
Part No.: 054101

1. 3-Octanone
2. *trans*-beta-Ocimene
3. Linalool L
4. Octenyl Acetate
5. 1-Borneol
6. Terpinene-4-ol
7. Linalyl Acetate
8. Lavandulyl Acetate
9. Santalene
10. *trans*-Caryophyllene


Note: The 5% phenyl BPX5 column is the perfect partner with a BP20 column for essential oil analysis.



ANALYSIS OF ROSEMARY ESSENTIAL OIL

Phase: BPX5, 0.25µm film
Column: 30 x 0.25mm ID
 Initial Temp.: 45°C, 1 min
 Rate 1: 5°C/min
 Final Temp.: 250°C, 15 min
 Detector: 5973 MSD
 Injector Mode: Split (50:1)
 Injector Temp.: 250°C
 Carrier Gas: Helium, 35 cm/sec @ 45°C
 Injection Volume: 1µL
 Concentration (ng/mL): 1:20 dilution in DCM
 GC model: HP 6890 GC, 5973 MSD
 Interface/Source/Quad: 320°C, 250°C, 106°C
Part No.: 054101

1. alpha-Pinene
2. Camphene
3. beta-Pinene
4. beta-Myrcene
5. p-Cymene
6. Limonene
7. 1,8-Cineole
8. gamma-Terpinene
9. Linalool L
10. Camphor
11. 1-Borneol
12. Terpinene-4-ol
13. alpha-Terpineol
14. *trans*-Caryophyllene




ANALYSIS OF EUCALYPTUS ESSENTIAL OIL

Phase: BPX5, 0.25µm film
Column: 30m x 0.25mm ID
 Initial Temp.: 45°C, 1 min
 Rate 1: 5°C/min
 Final Temp.: 250°C, 15 min
 Detector: 5973 MSD
 Injector Mode: Split (50:1)
 Injector Temp.: 250°C
 Carrier Gas: Helium, 35 cm/sec @ 45°C
 Injection Volume: 1mL
 Concentration (ng/mL): 1:20 dilution in DCM
 GC model: HP 6890 GC, 5973 MSD
 Interface/Source/Quad: 320°C, 250°C, 106°C
Part No.: 054101

1. alpha-Pinene
2. Limonene
3. 1,8-Cineole
4. Terpinene-4-ol
5. alpha-Terpineol
6. Camphene

Note: The non-polar nature of this phase will generally give separation based on boiling point. This is complimentary to the unique selectivity from the polyethylene glycol BP20 phase.




PHARMACEUTICALS - BASIC DRUG SCREEN (10-20ng/component)

Phase: BPX5, 1.0µm
Column: 25m x 0.53mm ID
 Initial Temp.: 120°C
 Rate: 10°C/min
 Final Temp.: 310°C
 Detector: FID
 Injector: Split, 240°C
 Carrier Gas: H₂, 2psi
Part No.: 054131

1. Methamphetamine
2. Phendimetrazine
3. Phencyclidine
4. Mepivocaine
5. Methaqualone
6. Amitriptyline
7. Codeine
8. Diazepam
9. Fentanyl


Note: A high degree of inertness and the high maximum operating temperature (for semi-volatile drugs), makes the BPX5 column suitable for this class of compounds.



STERIOD ANALYSIS (UNDERIVATIZED)

Phase: BPX5, 0.25µm
Column: 25m x 0.22mm ID
 Initial Temp.: 180°C
 Rate: 8°C/min
 Final Temp.: 350°C, 10 min
 Detector: FID
 Sensitivity: 32 x 10⁻¹² AFS
 Injection Mode: Split
 Carrier Gas: H₂, 10psi
Part No.: 054113

1. Testosterone
2. Pregnenolone
3. Progesterone
4. Cholesterol




Note: Whether your particular application is forensic screening or pharmaceutical testing, the BPX5 phase is ideal.

PHTHALATES - ANALYSIS OF PHTHALATE ESTERS

Phase: BPX5, 0.5µm film
Column: 25m x 0.32mm ID
 Initial Temp.: 180°C, 0min
 Rate: 15°C/min
 Final Temp.: 350°C, 5min
 Carrier Gas: H₂, 6psi
 Detector: FID
 Sensitivity: 140 x 10⁻¹² AFS
 Injection Mode: Split
Part No.: 054125


1. Dimethyl Phthalate
2. Diethyl Phthalate
3. Diisopropyl Phthalate
4. Dipropyl Phthalate
5. Diisobutyl Phthalate
6. Dibutyl Phthalate
7. Dipropyl Phthalate
8. Dihexyl Phthalate
9. Diheptyl Phthalate
10. Branched Dinonyl Phthalate
11. Dioctyl Phthalate
12. Diundecyl Phthalate
13. Didodecyl Phthalate



PETROLEUM - LEADED GASOLINE

Phase: BPX5, 0.25µm film
Column: 25m x 0.22mm ID
 Initial Temp.: -20°C, 1min
 Rate 1: 4°C/min
 Temp. 2: 40°C
 Rate 2: 8°C/min
 Final Temp.: 150°C
 Detector: FID
 Sensitivity: 10 x 10⁻¹² AFS
 Injection Mode: Split
 Carrier Gas: H₂, 10psi
Part No.: 054113

- | | | |
|-----------------------|----------------------------|-----------------------------|
| 1. iso-Butane | 12. 2,3-Dimethylpentane | 23. 1-Methyl-3-ethylbenzene |
| 2. n-Butane | 13. 3-Methylhexane | 24. 1-Methyl-4-ethylbenzene |
| 3. iso-Pentane | 14. 2,2,4-Trimethylpentane | 25. 1,3,5-Trimethylbenzene |
| 4. n-Pentane | 15. n-Heptane | 26. 1-Methyl-2-ethylbenzene |
| 5. 2,2-Dimethylbutane | 16. Toluene | 27. 1,2,4-Trimethylbenzene |
| 6. 2,3-Dimethylbutane | 17. 2,3-Dimethylhexane | 28. 1,2,3-Trimethylbenzene |
| 7. 2-Methylpentane | 18. Ethylbenzene | 29. Naphthalene |
| 8. 3-Methylpentane | 19. m-Xylene | 30. 2-Methylnaphthalene |
| 9. n-Hexane | 20. p-Xylene | 31. 1-Methylnaphthalene |
| 10. Benzene | 21. o-Xylene | |
| 11. 2-Methylhexane | 22. n-Propylbenzene | |



ORDERING INFORMATION - BPX5

Column Type	ID (mm)	Length (m)	Film Thickness (µm)	Temperature Limits (°C)	Part No.	
BPX5	0.10	10	0.10	-40 to 360/370	054099	
		10	1.20	-40 to 360/370	054106	
	0.15	12	0.25	-40 to 360/370	054103	
		12	0.40	-40 to 360/370	054107	
		25	0.25	-40 to 360/370	054104	
		25	0.40	-40 to 360/370	054108	
		30	0.15	-40 to 360/370	054110	
		50	0.25	-40 to 360/370	054105	
		0.22	12	0.25	-40 to 360/370	054112
			25	0.25	-40 to 360/370	054113
			25	1.00	-40 to 360/370	054116
			30	0.25	-40 to 360/370	054142
	50		0.25	-40 to 360/370	054114	
	0.25	50	1.00	-40 to 360/370	054117	
		15	0.25	-40 to 360/370	054100	
		15	1.00	-40 to 360/370	054121	
		30	0.25	-40 to 360/370	054101	
		30	0.50	-40 to 360/370	0541025	
		30	1.00	-40 to 360/370	054122	
		60	0.25	-40 to 360/370	054102	
	0.32	60	1.00	-40 to 360/370	054123	
		6	1.00	-40 to 360/370	0541261	
	0.53	12	0.25	-40 to 360/370	054118	
		12	0.50	-40 to 360/370	054124	
		12	1.00	-40 to 360/370	054127	
		15	0.25	-40 to 360/370	054144	
		15	1.00	-40 to 360/370	054152	
		25	0.25	-40 to 360/370	054119	
		25	0.50	-40 to 360/370	054125	
		25	1.00	-40 to 360/370	054128	
		25	3.00	-40 to 350/360	054136	
		30	0.25	-40 to 360/370	054145	
		30	0.50	-40 to 360/370	0541205	
		30	1.00	-40 to 360/370	054153	
		50	0.25	-40 to 360/370	054120	
		50	0.50	-40 to 360/370	054126	
		50	1.00	-40 to 360/370	054129	
		60	0.25	-40 to 360/370	054146	
		60	1.00	-40 to 360/370	054154	
		12	0.25	-40 to 360/370	054133	
		12	1.00	-40 to 360/370	054130	
		12	3.00	-40 to 350/360	054138	
		15	0.50	-40 to 360/370	0541344	
		15	1.00	-40 to 360/370	054147	
		15	1.50	-40 to 350/360	0541347	
	15	3.00	-40 to 350/360	054159		
	25	0.25	-40 to 360/370	054134		
	25	1.00	-40 to 360/370	054131		
	25	3.00	-40 to 350/360	054139		
	30	0.50	-40 to 360/370	0541345		
	30	1.00	-40 to 360/370	054148		
	30	1.50	-40 to 350/360	0541348		
	30	3.00	-40 to 350/360	054160		
50	1.00	-40 to 360/370	054132			
60	1.00	-40 to 360/370	054158			



www.sge.com

*Mentioned trademarks are registered trademarks of their respective owners.
Fused Silica capillary columns are manufactured under license granted by Agilent Technologies.

SGE International Pty Ltd
Toll Free: 1800 800 167
Tel: +61 (0) 3 9837 4200
Fax: +61 (0) 3 9874 5672
email: support@sge.com

SGE Europe Ltd (UK)
Tel: +44 (0) 1908 568844
Fax: +44 (0) 1908 566790
email: uk@sge.com

SGE (Italia) srl
Tel: +39 (06) 4429 0206
Fax: +39 (06) 4429 0724
email: sge.italia@tin.it

China Service Centre
Tel: +86 (10) 6588 8666
Fax: +86 (10) 6588 6577

SGE Laboratory Accessories Pvt Ltd (India)
Tel: +91 (022) 471 5896
Fax: +91 (022) 471 6592
email: sgeindia@vsnl.com

SGE Incorporated (USA)
Toll Free: 800 945 6154
Tel: +1 512 837 7190
Fax: +1 512 836 9159
email: usa@sge.com

SGE Europe Ltd (France)
Tel: +33 (0) 1 69 29 80 90
Fax: +33 (0) 1 69 29 09 25
email: france@sge.com

SGE Japan Inc
Tel: +81 (045) 222 2885
Fax: +81 (045) 222 2887
email: japan@sge.com

SGE GmbH (Germany)
Tel: +49 (0) 6151 860486
Fax: +49 (0) 6151 860489
email: germany@sge.com

©Copyright SGE International Pty Ltd. 2004. All rights reserved.
PD-0214-C Rev:01 07/04